3D PrintedCoCr Crowns & Bridges

01.

Highly accurate

Accurate fit with accuracy error of less than 0.02mm.





02.

Great bond strength

Outstanding bond strength exceeding 37MPa, ensuring a robust connection with your dental work.

03.

Minimal ion release & good biocompatibility

Minimized metal ion release, with levels consistently below 0.8ug/cm, realizing highest level of biocompatibility and patient safety.





Doubled Design Efficiency:

Metal printing improves design efficiency, enabling an average of 100 dental crowns to be designed by a single person in just 8 hours, doubling productivity.

03.

Enhanced Polishing Efficiency:

Metal 3D printing technology enhances polishing efficiency. On average, a single person can polish 70 dental crowns in 8 hours, which is 50% improvement in productivity.

02.

Faster Production:

Metal 3D printing and post-processing treatments, such as annealing, allow for the production of 100 dental crowns in just 7 hours. This is achieved without the need for complex and timeconsuming steps like re-moldling, embedding, was loss, or casting.

Superior Quality Control:

The quality control achieved through metal 3D printing is exceptional, with a 99.9% yield rate of defect-free products. Dental crowns produced through this method exhibit a flawless surface with no air bubbles or sand inclusions. They remain structurally stable with no deformation, and there are no casting defects or missing components.







